

GABARAP

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 13kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GABARAP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ch

Synonym:

MM46; ATG8A; GABARAP-a;

Catalog #: A5616

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11337

Isotype: IgG

Swiss Prot: O95166

Purity: Affinity purification

For research use only.

Background:

Gamma-aminobutyric acid A receptors

To place an order, please [Click HERE](#).