

MYL6

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MYL6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LC17; ESMLC; LC17A; LC17B; MLC-3; MLC1SM; MLC3NM; MLC3SM; LC17-GI; LC17-NM

Catalog #: A8557

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4637

Isotype: IgG

Swiss Prot: P60660

Purity: Affinity purification

For research use only.

Background:

Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain that is expressed in smooth muscle and non-muscle tissues. Genomic sequences representing several pseudogenes have been described and two transcript variants encoding different isoforms have been identified for this gene.

To place an order, please [Click HERE](#).