

## Der P1

**Description:** The E.Coli derived recombinant protein contains the Dermatophagoides pteronyssinus Der P1 full length protein (a.a. 1-321) and fused to a 6 His Tag at C-terminus, having a total Mw of 34.5kDa, pI 5.6.

**Catalog #:**ALPS-010

**Synonyms:**Peptidase 1, Major mite fecal allergen Der p 1, Allergen Der p I, Der p 1, DERP1, Der-P1.

For research use only.

**Purity:**Protein is >95% pure as determined by 10% SDS-PAGE (coomassie staining).

**Formulation:**

(1mg/ml) 60mM NaCl, 50mM Tris-HCl pH 8.0 and 1.2M Urea.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

DERP1 is a thiol protease, with a preference for substrates with a large hydrophobic side chain in the P2 position, or with basic residues. DERP1 is a C1 peptidase family member. DERP1 has extensive endopeptidase specificity. DERP1 is N-glycosylated. N-glycanase treatment does not completely remove carbohydrates, suggesting that the protein contains additional glycosylation sites. DERP1 causes an allergic reaction in humans. Common symptoms of mite allergy are bronchial asthma, allergic rhinitis and conjunctivitis. DERP1 binds to IgE in 80% of patients with house dust allergy.

**Storage:**

Der-P1 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

**To place an order, please [Click HERE](#).**