

## Phospho-Src-Y527

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 60kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding Y527 of human Src

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

oqt

**Synonym:**

ASV; SRC1; c-SRC; p60-Src;

**Catalog #:** AP0470

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 6714

**Isotype:** IgG

**Swiss Prot:** P12931

**Purity:** Affinity purification

For research use only.

**Background:**

This gene is highly similar to the v-src gene of Rous sarcoma virus. This proto-oncogene may play a role in the regulation of embryonic development and cell growth. The protein encoded by this gene is a tyrosine-protein kinase whose activity can be inhibited by phosphorylation by c-SRC kinase. Mutations in this gene could be involved in the malignant progression of colon cancer. Two transcript variants encoding the same protein have been found for this gene.

*To place an order, please [Click HERE](#).*