

CXCL4 Variant 1 Human

Description: CXCL4 Variant-1 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 77 amino acids and having a molecular mass of 8.7 kDa. The CXCL4 Variant-1 is fused to 6xHis tag at N-Terminus and purified by standard chromatography techniques.

Catalog #: CHPS-250

For research use only.

Synonyms: CXCL4, PF-4, PF4, Iroplact, Oncostatin-A, SCYB4, MGC138298.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered white lyophilized powder.

Amino Acid Sequence: MHHHHHHEAE EDGDLQCLCV KTTSQVRPRH ITSLEVIKAG
PHCPTAQLIA TLKNGRKICL DLQALLYKKI IKEHLES.

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

The protein was lyophilized without additives.

Stability:

Human CXCL4 although stable at 25°C 1 week, should be stored desiccated below -18°C. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized CXCL4 in sterile 18M-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Introduction:

Platelet factor-4 is a 70-amino acid protein that is released from the alpha-granules of activated platelets and binds with high affinity to heparin. Its major physiologic role appears to be neutralization of heparin-like molecules on the endothelial surface of blood vessels, thereby inhibiting local antithrombin III activity and promoting coagulation. As a strong chemoattractant for neutrophils and fibroblasts, PF4 probably has a role in inflammation and wound repair. Oncostatin-A is a member of the CXC chemocine family. Human PF4 is used for the proof of heparin-induced thrombocytopenia. Furthermore it is used as an inhibitor in the angiogenesis during tumor therapy.

To place an order, please [Click HERE](#).