

## IFNW1 Human

**Description:** Interferon-Omega 1 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 172 amino acids and having a molecular mass of 19.9kDa. The Interferon-Omega 1 is purified by proprietary chromatographic techniques.

**Catalog #:** CYPs-047

For research use only.

**Synonyms:** Interferon omega-1, Interferon alpha-II-1, IFNW1.

**Source:** Escherichia Coli.

**Physical Appearance:** Sterile Filtered White lyophilized (freeze-dried) powder.

**Amino Acid Sequence:** MCDLPQNHGL LSRNTLVLLH QMRRISPFLC LKDRRDFRFP  
QEMVKGSQLQ KAHVMSVLHE MLQQIFSLFH TERSSAAWNM TLLDQLHTGL HQQLQHLETC  
LLQVVGEGES AGAISSPALT LRRYFQGIRV YLKEKKYSDC AWEVVRMEIM KSLFLSTNMQ  
ERLRSKDRDL GSS.

**Purity:** Greater than 95.0% as determined by: (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

**Formulation:**

Lyophilized from a 0.2

**Stability:**

Lyophilized IFNW1 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution IFNW1 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Solubility:**

It is recommended to reconstitute the lyophilized IFNW1 in sterile 18M-cm H<sub>2</sub>O not less than 100

**Introduction:**

IFN-Omega 1 is a type I interferon, that can be induced by virus-infected leukocytes. Type I interferon family members, which include IFN-alpha, IFN-beta, and IFN-omega, signal through IFNAR-1/IFNAR-2 receptor complex, and wield antiviral and antiproliferative activities. IFNW1 exhibits about 75% sequence homology with IFN-a, and contains 2 conserved disulfide bonds, which are essential for full biological activity.

**Biological Activity:**

The ED<sub>50</sub> was determined by a cytotoxicity assay using Human TF-1 cells is < 0.01 ng/ml, corresponding to a specific activity of > 100,000,000 units/mg.

**To place an order, please [Click HERE](#).**