

## Leptin Receptor Chicken

**Description:**Leptin Binding Domain Chicken Recombinant also called Leptin Receptor produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 208 amino acids and having a molecular mass of 24.5 kDa. Chicken Leptin Receptor consists of the cytokine binding domain of leptin receptor amino acids 420-626 of chicken leptin receptor. The Leptin Binding Domain is purified by proprietary chromatographic techniques.

**Catalog #:**CYPS-516

For research use only.

**Synonyms:**OB Protein, Obesity Protein, OBS, Obesity factor, Leptin Receptor.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile Filtered colorless solution at a concentration of 0.4 mg/ml.

**Amino Acid Sequence:**The sequence of the first six N-terminal amino acids was determined and was found to be Ala-Ile-Asp-Val-Asn-Ile. Biological Activity is evidenced by high affinity binding of mammalian leptins at 1:1 molar ratio.

**Purity:**Greater than 99.0% as determined by: (a) Analysis by SEC-HPLC. (b) Analysis by SDS-PAGE.

**Formulation:**

The protein was filter sterilized and stored at 4°C (0.2 to 0.5 mg/ml) solution of Tris-HCl buffer, pH 9.0 with 150mM NaCl.

**Stability:**

Sterile solutions at 0.5mg/ml or less are stable at 4°C for several months.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**To place an order, please [Click HERE](#).**