

## NUDT16L1 Human

**Description:**NUDT16L1 Human Recombinant fused with a 20 amino acid His tag at N-terminus produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 231 amino acids (1-211 a.a.) and having a molecular mass of 25.5kDa. The NUDT16L1 is purified by proprietary chromatographic techniques.

**Catalog #:**ENPS-086

For research use only.

**Synonyms:**Protein syndesmos, NUDT16-like protein 1, NUDT16L1, SDOS, MGC11275.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile Filtered colorless solution.

**Amino Acid Sequence:**MGSSHHHHHH SSGLVPRGSH MSTAAVPELK QISRVEAMRL  
GPGWSHSCHA MLYAANPGQL FGRIPMRFSV LMQMRFDGLL GFPGGFVDRR FWSLEDGLNR  
VLGLGLGCLR LTEADYLSSH LTEGPHRVVA HLYARQLTLE QLHAVEISAV HSRDHGLEVL  
GLVRVPLYTQ KDRVGGFPNF LSNAFVSTAK CQLLFALKVL NMMPEEKLVE ALAAATEKQK  
KALEKLLPAS S.

**Purity:**Greater than 95.0% as determined by SDS-PAGE.

**Formulation:**

The NUDT16L1 solution (0.5 mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 1mM DTT.

**Stability:**

NUDT16L1 should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

Syndesmos (NUDT16L1) is a cytoplasmic protein which interacts specifically with the cytoplasmic domain of syndecan-4, and co-localizes with syndecan-4 in focal contacts. NUDT16L1 interacts with the focal adhesion adaptor protein paxillin. The binding of NUDT16L1 to paxillin is direct, and these interactions are triggered by the activation of protein kinase C.

**To place an order, please [Click HERE](#).**