

## BCOADC-E2 Human

**Description:**BCOADC-E2 is a Full-length cDNA coding for the mature form of the human BCOADC-E2 protein having a molecular mass of 47,321 Dalton(pH 6.5). BCOADC-E2 protein is fused to a hexa-histidine purification tag.

**Catalog #:**ENPS-091

**Source:**Sf9 insect cells.

For research use only.

**Purity:**Greater than 75% as determined by SDS-PAGE.

**Formulation:**

DLST (0.91mg/ml) is supplied in 20mM HEPES buffer pH-8.0, 200mM NaCl and 20% glycerol.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Applications:**

Western blot with anti-M2-Antigen autoantibody-positive patient sera or monoclonal anti-hexa-His-tag antibody.

**Storage:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

**To place an order, please [Click HERE](#).**