

MTHFS Human

Description: MTHFS produced in E.Coli is a single, non-glycosylated polypeptide chain containing 223 amino acids (1-203a.a.) and having a molecular mass of 25.4 kDa. MTHFS is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: ENPS-103

For research use only.

Synonyms: 5,10-methenyltetrahydrofolate synthetase (5-formyltetrahydrofolate cyclo-ligase), HsT19268, Methenyl-THF synthetase, FLJ30410, EC 6.3.3.2.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MAAAVSSAK RSLRGELKQR
LRAMSAEERL RQSRVLSQKV IAHSEYQKSK RISIFLSMQD EIETEEIIKD IFQRGKICFI
PRYRFQSNHM DMVRIESPEE ISLLPKTSWN IPQPGEGDVR EEALSTGGLD LIFMPGLGFD
KHGNRLGRGK GYYDAYLKRC LQHQEVKPYT LALAFKEQIC LQVPVNENDM KVDEVLYEDS
STA

Purity: Greater than 95% as determined by SDS-PAGE.

Formulation:

The MTHFS protein solution (0.5mg/1ml) is formulated in 20mM Tris-HCl buffer (pH8.0), 200mM NaCl, 5mM DTT and 30% glycerol.

Usage:

NeoBiolabs products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

MTHFS is a cytosolic protein which takes part in the formate metabolic process. MTHFS along with a magnesium cofactor catalyzes the ATP-dependent reaction which reduces 5-formyltetrahydrofolate (5-MTHF) to 5,10-methenyltetrahydrofolate (MTHF). MTHF is the substrate used by MTHFR (methylenetetrahydrofolate reductase) to generate 5-MTHF. In addition, MTHF is a coenzyme used in thymidine biosynthesis by thymidylate synthase (FAD).

Storage:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

To place an order, please [Click HERE](#).