

HDAC2 Human

Description: HDAC2 Human Recombinant produced in Hi-5 Cell is a single, non-glycosylated polypeptide chain containing 496 amino acids (1-488) and having a molecular mass of 56.4 kDa. The HDAC2 is fused to an 8 amino acid His-Tag at C-terminus and purified by proprietary chromatographic techniques.

Catalog #: ENPS-164

For research use only.

Synonyms: Histone deacetylase 2, YAF1, HD2, YY1-associated factor 1, transcriptional regulator homolog RPD3, RPD3, EC 3.5.1.98.

Source: Hi-5 Cell.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MAYSQGGGKK KVCYYYDGD I GNYYYGQGHP MKPHRIRMTH
NLLLNLYGLYR KMEIYRPHKA TAEEMTKYHS DEYIKFLRSI RPDNMSEYSK QMQRFNVGED
CPVFDGLFEF CQLSTGGSVA GAVKLNRRQT DMAVNWAGGL HHAKKSEASG FCYVNDIVLA
ILELLKYHQR VLYIDIDIHHGDGVVEAFYT TDRVMTVSFH KYGEYFPGTG DLRDIGAGKG
KYYAVNFPMR DGI

Purity: Greater than 85% as determined by SDS-PAGE.

Formulation:

The HDAC2 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 0.1M NaCl, 0.1mM PMSF and 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HDAC2 is a member of the histone deacetylase family that performs through the construction of large multiprotein complexes and are in charge of the deacetylation of lysine residues on the N-terminal region of the core histones. HDAC2 forms transcriptional repressor complexes by relating to a diversity of proteins, like YY1- a mammalian zinc-finger transcription factor. In addition, HDAC2 has a vital part in transcriptional regulation, cell cycle progression and developmental procedures.

To place an order, please [Click HERE](#).