

## PSMB4 Human

**Description:** PSMB4 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 240 amino acids (46-264 a.a.) and having a molecular mass of 26.6kDa. PSMB4 is fused to a 21 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

**Catalog #:** ENPS-204

For research use only.

**Synonyms:** Proteasome (prosome, macropain) subunit beta type 4, PROS26, HN3, Macropain beta chain, Multicatalytic endopeptidase complex beta chain, Proteasome chain 3, 26 kDa prosomal protein, HsN3, PROS-26, hsBPROS26, EC 3.4.25.1.

**Source:** Escherichia Coli.

**Physical Appearance:** Sterile Filtered colorless solution.

**Amino Acid Sequence:** MGSSHHHHHH SSGLVPRGSH MTQNPMTGT SVLGVKFEGG  
VVIAADMLGS YGSLARFRNI SRIMRVNNST MLGASGDYAD FQYLKQVLGQ MVIDEELLGD  
GHSYSPIAIH SWLTRAMYSR RSKMNPLWNT MVIGGYADGE SFLGYVDMLG VAYEAPSLAT  
GYGAYLAQPL LREVLEKQPV LSQTEARDLV ERCMRVLYR DARSYNRFQI ATVTEKGVEI  
EGPLSTETNW DI

**Purity:** Greater than 90.0% as determined by SDS-PAGE.

**Formulation:**

PSMB4 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 100mM NaCl, 1mM DTT and 30% glycerol.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

PSMB4 is a member of the peptidase T1B family. The 20S Proteasome chamber holds structural alpha subunits and predominantly catalytic beta subunits. The two external rings in the proteasome contain seven alpha subunits each and the two internal rings contain of seven beta subunits each. PSMB4 is a beta subunit of the 20S Proteasome that is spread all over the eukaryotic cells in large quantities and cleaves peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway.

**To place an order, please [Click HERE](#).**