

## PSMA7 Human

**Description:** PSMA7 Human Recombinant fused to N-terminal His-Tag produced in E.Coli is a single, non-glycosylated polypeptide chain containing 268 amino acids (1-248) and having a molecular mass of 30kDa. PSMA7 is fused to a 20 amino acid His Tag at N-Terminus and purified by standard chromatography techniques.

**Catalog #:** ENPS-410

For research use only.

**Synonyms:** HSPC, RC6-1, XAPC7, MGC3755, PSMA-7, Proteasome subunit alpha type-7, Proteasome subunit RC6-1, Proteasome subunit XAPC7, PSMA7, C6.

**Source:** Escherichia Coli.

**Physical Appearance:** Sterile Filtered colorless solution.

**Amino Acid Sequence:** MGSSHHHHHH SSGLVPRGSH MSYDRAITVF SPDGHILFQVE  
YAEAVKKS TAVGVRGRDI VVLGVEKKS AKLQDERTVR KICALDDNVC MAFAGLTADA  
RIVINRARVE CQSHRLTVED PVTVEYITRY IASLKQRYTQ SNGRRPFGIS ALIVGFDFDG  
TPRLYQTDPS GTYHAWKANA IGRGAKSVRE FLEKNYTDEA IETDDLTIKL VIKALLEVVQ  
SGGKNIELAV MR

**Purity:** Greater than 90.0% as determined SDS-PAGE.

**Formulation:**

The PSMA7 solution contains 20mM Tris-HCl pH-8, 0.1M NaCl, and 20% glycerol.

**Stability:**

PSMA7 although stable 4°C for 4 weeks, should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

The proteasome is a multicatalytic proteinase complex with a highly assembled ring-shaped 20S core structure which is composed of 4 rings of 28 non-identical subunits, 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. Proteasomes are found throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. PSMA7 is part of the peptidase T1A family, that is a 20S core alpha subunit. PSMA7 interacts particularly with the hepatitis B virus X protein, a protein critical to viral replication. PSMA7 is involved in regulating hepatitis virus C internal ribosome entry site activity that is crucial for viral replication. PSMA7 is in charge for regulating the hypoxia-inducible factor-1alpha, a transcription factor important for cellular responses to oxygen tension. PSMA7 is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH. PSMA7 interacts specifically with two subdomains of HIF-1alpha and inhibited the transactivation function of HIF-1alpha under both normoxic and hypoxia-mimicking conditions.

**To place an order, please [Click HERE](#).**