

Thymosin 4

Description:Thymosin b4 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 43 amino acids and having a molecular mass of 4.9kDa.The Thymosin b4 Human is purified by proprietary chromatographic techniques.

Catalog #:HOPS-010

For research use only.

Synonyms:Thymosin beta-4, T beta 4, Fx , TB4X, PTMB4, TMSB4.

Source:EscherichiaColi.

Physical Appearance:SterileFiltered White lyophilized (freeze-dried) powder.

Amino Acid Sequence:SDKPDMAEIEKFDKSKLKKT ETQEKNPLPS KETIEQEKQA GES

Purity:Greater than 95.0% as determined by SDS-PAGE and HPLC analyses.

Formulation:

Lyophilized from a 0.2um filtered concentrated solution in 1xPBS, pH 7.4.

Stability:

Lyophilized Thymosin B4 Human although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Thymosin B4 Human should be stored at 4°C between 2-7 days and for future use below -18°C. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized Thymosin B4 in sterile 18M-cm H2O not less than 100

Introduction:

CB1 cannabinoid receptor-interacting protein 1 (CNRIP1) is a 164 amino acid protein and G-protein coupled receptor which is a member of the CNRIP family. The CNRIP1 interacts with the C-terminal tail of cannabinoid receptor 1. CNRIP1 is expressed in the brain tissue and, at low levels, in the testis. CNRIP1 is involved in appetite, synaptic plasticity, neuroprotection and analgesia.

To place an order, please [Click HERE](#).