

LH Porcine

Description: Porcine Luteinizing Hormone is a glycosylated, polypeptide chain which stimulates maturation of follicle, induces ovulation, accelerates formation of corpus luteum and secretion of pregnanolone. Porcine Luteinizing Hormone delays ovulation and ovarian cyst.

Catalog #:HOPS-317

Source: Porcine.

For research use only.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation:

The LH was lyophilized with no additives.

Stability:

Lyophilized Porcine LH although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution LH should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized Porcine LH in sterile 18M-cm H₂O not less than 100

To place an order, please [Click HERE](#).