

HBsAg ayw

Catalog #:HXPS-877

Source:Saccharomyces cerevisiae.

Physical Appearance:Sterile Filtered clear solution.

For research use only.

Purity:Greater than 98.0% as determined by SDS-PAGE.

Formulation:

Sterile Filtered solution containing 0.05M phosphate buffer pH 7.2 and 0.2M NaCl.

Stability:

HBsAg should be stored at 4°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV). The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.Recombinant HbsAg ayw full length is a 24kDa protein cloned from HBV 320 genome.

To place an order, please [Click HERE](#).