

HBsAg adr CHO

Description: Recombinant HBsAg adr produced in CHO cells, the molecular weight is approximately 23-27kDa.

Catalog #: HXPS-882

Source: CHO.

For research use only.

Physical Appearance: Sterile Filtered clear solution.

Purity: Greater than 95.0% as determined by RP-HPLC.

Formulation:

Sterile Filtered solution containing 20mM PB and 154mM NaCl.

Stability:

HBsAg Should be stored at 4°C, DO NOT FREEZE.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV) S-gene. The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.

To place an order, please [Click HERE](#).