

HBsAg ayw

Description: Recombinant HbsAg ayw full length is a 31kDa protein cloned from HBV 320 genome and fused to 6xHis Tag.

Catalog #: HXPS-884

Source: Escherichia Coli.

For research use only.

Physical Appearance: Sterile Filtered clear solution.

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

HBsAg ayw solution containing 7.5 mM phosphate buffer pH-7.2, 75 mM NaCl and 50% glycerol.

Stability:

HBsAg Should be stored at 4°C between 3-4 weeks and for future use should be stored below -18°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV). The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.

To place an order, please [Click HERE](#).