

MV Hemagglutinin Mosaic

Description: The E.Coli derived recombinant protein contains the large hemagglutinin immunodominant regions 1-30, 115-150 and 379-410 amino acids.

Catalog #:MEPS-012

Purity: Measles Mosaic protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

For research use only.

Purification Method:

Measles Mosaic protein was purified by proprietary chromatographic technique.

Formulation:

50mM Tris-HCl, pH 8.0, 60mM NaCl, 10mM glutathione and 50% glycerol.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Measles Mosaic protein can be used as an antigen in ELISA and Western Blots. Immunoreactive with sera of Measles Virus infected individuals.

Introduction:

Measles, also known as rubeola, caused by a virus, specifically a paramyxovirus of the genus Morbillivirus. To date, 21 strains of the measles virus have been identified. Measles virus (MV) is an enveloped, nonsegmented negative-stranded RNA virus of the Paramyxoviridae family.

Storage:

Measles Mosaic although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

To place an order, please [Click HERE](#).