

CSNK2A1 Human, His

Description: CK2A Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain encoding the sequence of 411 amino acids and having a molecular mass of 47.3 kDa. Casein Kinase 2 alpha subunit is purified by proprietary chromatographic techniques.

Catalog #: PKPS-232

Synonyms: Casein kinase II subunit alpha, EC 2.7.11.1, CK II, CK2A1, CKII alpha, CSNK2A1, PKCK2, CSK21.

For research use only.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MSGPVPSRAR VYTDVNTHRP
REYWDYESHV VEWGNQDDYQLVRKLGRGKY SEVFEAINIT NNEKVVKIL KPVKKKKIKR
EIKILENLRG GPNIITLADI VKDPVSRTPA LVFEHVNNTD FKQLYQTLTD YDIRFYMYEI
LKALDYCHSM GIMHRDVKPH NVIMIDHEHRK LRLIDWGLAE FYHPGQEYNV RVASRYFKGP
ELLVDYQMYD YSL

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

CSNK2A1 is supplied in 20mM Tris pH-8.0 , 500mM NaCl, 1mM DTT and 50% glycerol.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

The Casein kinase 2 (EC2.7.11.1) is a serine/threonine-selective protein kinase that is a tetramer of two alpha subunits and two beta subunits. The alpha subunits have the catalytic kinase domain. Casein kinase 2 has been implicated in cell cycle control, DNA repair, regulation of the circadian rhythm and other cellular processes. Casein kinase 2 activity has been reported to be activated following Wnt signaling pathway activation. A Pertussis toxin-sensitive G protein and Dishevelled appear to be an intermediary between Wnt-mediated activation of the Frizzled receptor and activation of casein kinase 2. Mice that lack casein kinase 2 alpha prime have a defect in the morphology of developing sperm.

Storage:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

To place an order, please [Click HERE](#).