

RAB5B Human

Description: RAB5B Human Recombinant fused with a 23 amino acid His tag at N-terminus produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 238 amino acids (1-215 a.a.) and having a molecular mass of 26.1kDa. The RAB5B is purified by proprietary chromatographic techniques.

Catalog #: PRPS-063

For research use only.

Synonyms: Ras-related protein Rab-5B, RAB5B.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGSMTSRSTA RPNGQPQASK
ICQFKLVLLG ESAVGKSSLV LRFVKGQFHE YQESTIGAAF LTQSVCLDDT TVKFEIWDTA
GQERYHSLAP MYYRGAQAAI VVYDITNQET FARAKTWVKE LQRQASPSIV IALAGNKADL
ANKRMVEYEE AQAYADDNSL LFMETSAKTA MNVNDLFLAI AKKLPKSEPQ NLGGAAGRSR
GVDLHEQSQQ NK

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

The RAB5B solution (1 mg/ml), 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 20% glycerol and 0.2M NaCl.

Stability:

RAB5B should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Ras-related protein Rab-5B (RAB5B) belongs to the Rab family of small (monomeric) G proteins. RAB5B is involved in endocytosis and recycling of cell surface molecules. RAB5B interacts with RIN2 and RIN3, which regulate its function, possibly by acting as GEFs. Knockdown of RAB5B abolished group I metabotropic glutamate receptor (mGluR)-mediated neuroprotection.

To place an order, please [Click HERE](#).