

C1QTNF8 Human

Description: C1QTNF8 Protein is a 28.4 kDa protein containing 255 aa fused to a 10 aa N-Terminal His-tag.

Catalog #: PRPS-141

Synonyms: CTRP8, UNQ5829, C1q/TNF-related protein 8, Complement C1q tumor necrosis factor-related protein 8.

For research use only.

Source: E. coli

Amino Acid Sequence: MKHHHHHHAS WPGLPRRPCV HCCRPWWPPG PYARVSDRDL
WRGDLWRGLP RVRPTIDIEI LKGEKGEAGV RGRAGRSGKE GPPGARGLQG RRGQKGQVGP
PGAACRRAYA AFSVGRRAYA AFSVGRREGL HSSDHFQAVP FDTLVNLDG AFDLAAGRFL
CTVPGVYFLS LNVHTWNYKE TYLHIMLNRR PAAVLYAQPS ERSVMQAQSL MLLLAAGDAV
WVRMFQRDRD NA

Formulation:

C1QTNF8 Human was filtered (0.4

Stability:

Store lyophilized C1QTNF8 Human at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles. Reconstituted C1QTNF8 can be stored at 4°C for a limited period of time; it does not show any change after two weeks at 4°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Western blotting

Solubility:

Add 0.1M Acetate buffer pH4 to prepare a working stock solution of approximately 0.5 mg/ml and let the lyophilized pellet dissolve completely. For conversion into higher pH value, we recommend intensive dilution by relevant buffer to a concentration of 10g/ml. In higher concentrations the solubility of this antigen is limited. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

To place an order, please [Click HERE](#).