

C16ORF53 Human

Description: C16ORF53 produced in E.Coli is a single, non-glycosylated polypeptide chain containing 274 amino acids (1-254 a.a.) and having a molecular mass of 29.9kDa (Molecular weight on SDS-PAGE will appear higher). C16ORF53 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-232

For research use only.

Synonyms: PAXIP1-associated protein 1, PTIP-associated protein 1, PA1, C16orf53, GAS.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MSLARGHGD T AASTAAPLSE
EGEVTSG LQA LAVEDTGGPS ASAGKADEG EGGREETERE GSGGEEAQGE VPSAGGEEPA
EEDSEDWCVP CSDEEVELPA DGQPWMPPPS EIQRLYELLA AHGTLELQAE ILPRRPPTPE
AQSEEERSDE EPEAKEEEEE KPHMPTEFDF DDEPVTPKDS LIDRRRTPGS SARSQKREAR
LDKVLSDMKR HK

Purity: Greater than 85.0% as determined by SDS-PAGE.

Formulation:

C16ORF53 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 10% glycerol, 2mM DTT and 0.1M NaCl.

Stability:

C16ORF53 Human Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

C16ORF53 is a component of a Set1-like multiprotein histone methyltransferase complex. C16ORF53 interacts with PAXIP1/PTIP; this interaction is direct and is necessary for the association with the rest of the PTIP complex. The C16ORF53 protein has a crucial role in maintaining genome stability, condensation of chromatin and progression through mitosis.

To place an order, please [Click HERE](#).