

CD34 Human, GST

Description: CD34 Human Recombinant full length expressed in E.coli, CD34 is antibody reactive. The CD34 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-299

Synonyms: Hematopoietic progenitor cell antigen CD34, CD34.

For research use only.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Formulation:

CD34 in 50mM Tris-HCl, pH7.5, 10mM L-glutathione (reduced).

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

CD34 is a single chain transmembrane glycoprotein which is selectively expressed on human lymphoid and myeloid hemapoeitic progenitor cells. CD34 has been used to measure angiogenesis, which reportedly predicts tumor recurrence.

To place an order, please [Click HERE](#).