

C1QTNF6 Human

Description: C1QTNF6 Human Recombinant produced in E.Coli is a single, non-glycosylated, Polypeptide chain containing 242 amino acids and having a molecular mass of 27 kDa. The protein contains an extra His tag at N-terminus. The C1QTNF6 amino acid sequence is identical to UniProtKB/Swiss-Prot entry Q9BXI9 amino acids 29259. The C1QTNF6 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-662

For research use only.

Synonyms: CTRP6, C1QTNF6, ZACRP6, Complement C1q tumor necrosis factor-related protein 6.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Amino Acid Sequence: MKHHHHHHAS TFDRAVASGC QRCCDSEDPL DPAHVSSASS
SGRPHALPEI RPYINITILK GDKGDPGPMG LPGYMGREGP QGEPGPQGSK GDKGEMGSPG
APCQKRFFAF SVGRKTALHS GEDFQTLLE RVFVNLDGCF DMATGQFAAP LRGIYFFSLN
VHSWNYKETY VHIMHNQKEA VILYAQPSE SIMQSQSVML DLAYGDRVVV RLFKRQRENA
IYSNDFDTYI TF

Purity: Greater than 95% as determined by SDS PAGE.

Formulation:

0.23mg/ml in 30mM Acetate Buffer pH-4.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

C1QTNF6 is a collagen type-IV specific for basement membranes. C1QTNF6 is vital for embryonic development.

To place an order, please [Click HERE](#).