

ATOX1 Human

Description: ATOX1 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 88 amino acids (1-68 a.a.) and having a molecular mass of 9.5kDa. ATOX1 is fused to a 20 amino acid His Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: PRPS-761

For research use only.

Synonyms: Antioxidant protein 1, ATX1, HAH1, Copper transport protein ATOX1, Metal transport protein ATX1, ATOX1, MGC138453, MGC138455.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MPKHEFSVDM TCGGCAEAVS
RVLNKLGGVK YDIDLPNKKV CIESEHSDMT LLATLKKTGK TVSYLGLE.

Purity: Greater than 95% as determined by SDS-PAGE.

Formulation:

ATOX1 solution containing 20mM Tris pH-8, 1mM DTT and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

ATOX1 is a metal transport protein that is part of the ATX1 family. ATOX1 is a copper chaperone that takes part in cellular antioxidant defense and can bind and deliver cytosolic copper to the copper ATPase proteins in the trans-Golgi network for later incorporation to the ceruloplasmin. ATOX1 plays a role as an antioxidant against superoxide and hydrogen peroxide, and consequently, takes an important part in cancer carcinogenesis. Because of ATOX1 cytogenetic location, the gene can be a good candidate gene for 5q-syndrome.

To place an order, please [Click HERE](#).