

## ANXA3 Human

**Description:**ANXA3 produced in E.Coli is a single, non-glycosylated polypeptide chain containing 323 amino acids (1-323 a.a.) and having a molecular mass of 36.3 kDa.ANXA3 is purified by proprietary chromatographic techniques.

**Catalog #:**PRPS-772

For research use only.

**Synonyms:**ANX3, ANXA3, Annexin-III, Lipocortin-3, Lipocortin-III, Annexin-3, Annexin A3, Placental anticoagulant protein III, PAP-III, 35-alpha calcimedlin, Inositol 1,2-cyclic phosphate 2-phosphohydrolase.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile filtered colorless solution.

**Amino Acid Sequence:**MASIWVGHRG TVRDYPDFSP SVDAEAIQKA IRGIGTDEKM  
LISILTERSN AQRQLIVKEY QAAYGKELKD DLKGDLSGHF EHLMVALVTP PAVFDAKQLK  
KSMKGAGTNE DALIEILTTR TSRQMKDISQ AYYTVYKKSL GDDISSETSG DFRKALLTLA  
DGRDES LKV DEHLAKQDAQ ILYKAGENRW GTDEDKFTEI LCLRSFPQLK LTFDEYRNIS  
QKDIVDSIKG EL

**Purity:**Greater than 90% as determined by SDS-PAGE.

### Formulation:

The ANXA3 protein solution contains 20mM Tris-HCl, pH-7.5, 1mM DTT and 10% glycerol.

### Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

### Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Introduction:

ANXA3 is an inhibitor of phospholipase A2. ANXA3 possesses anti-coagulant properties. ANXA3 cleaves the cyclic bond of inositol 1,2-cyclic phosphate to form inositol 1-phosphate. ANXA3 is part of the annexin family which are members of the calcium-dependent phospholipid-binding protein family that are involved in the regulation of cellular growth and in signal transduction pathways. ANXA3 takes part as an anti-coagulation factor.

**To place an order, please [Click HERE](#).**