

RAB11A Human

Description: RAB11A Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 249 amino acids(1-213a.a.) and having a molecular mass of 28.1 kDa. RAB11A protein is fused to a 36 amino acid His tag at N-terminus and is purified by standard chromatography.

Catalog #: PRPS-833

For research use only.

Synonyms: RAB11, YL8, Ras-related protein Rab-11A, RAB11A.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMGTR
DDEYDYLKFKV VLIGDSGVGK SNLLSRFTRN EFNLESKSTI GVEFATRSIQ VDGKTIKAQI
WDTAGQERYR AITSAYYRGA VGALLVYDIA KHLTYENVER WLKELRDHAD SNIVIMLVGN
KSDLRHLRAV PTDEARAF AE KNGLSFIETS ALDSTNVEAA FQTILTEIYR IVSQKQMSDR
RENDMSPSNN VV

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

RAB11A Human solution containing 20mM Tris-HCl pH-8, 1mM DTT, 1mM EDTA & 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

RAB11A is part of the small GTPase superfamily, Rab family. RAB11A is related with both constitutive and regulated secretory pathways, and takes part in protein transport. RAB11A is a well-known marker for protein trafficking, sorting, and recycling in the endosomal pathway. RAB11A is a main regulator of membrane delivery during cytokinesis.

To place an order, please [Click HERE](#).