

SARS Envelope

Description: The E.coli derived 34 kDa recombinant protein contains the N-terminus Enevelope protein 1-76 amino acids immunodominant regions.

Catalog #: SAPS-259

Purity: SARS Envelope protein is >95% pure as determined by 10% PAGE (coomassie staining).

For research use only.

Purification Method:

SARS Envelope protein was purified by proprietary chromatographic technique.

Specificity:

SARS Envelope protein is Immunoreactive with sera of SARS-infected individuals.

Formulation:

50mM Tris-HCl, 60mM NaCl, and 50% glycerol.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

SARS Envelope antigen is suitable for ELISA and Western blots, excellent antigen for detection of SARS with minimal specificity problems.

Introduction:

SARS Coronavirus is an enveloped virus containing three outer structural proteins, namely the membrane (M), envelope (E), and spike (S) proteins. Spike (S)-glycoprotein of the virus interacts with a cellular receptor and mediates membrane fusion to allow viral entry into susceptible target cells. Accordingly, S-protein plays an important role in virus infection cycle and is the primary target of neutralizing antibodies.

Storage:

SARS Envelope Protein is shipped at ambient temperature. Upon arrival, store at -20C.

To place an order, please [Click HERE](#).