

SARS Mosaic S(C)

Description: The E.Coli derived 37 kDa recombinant mosaic protein contains the C-terminal t section of the Spike protein 1051-1076, 1121-1154, 1162-1190 amino acids immunodominant regions.

Catalog #: SAPS-264

For research use only.

Purity: SARS Mosaic Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Purification Method:

SARS Mosaic protein was purified by proprietary chromatographic technique.

Specificity:

SARS Mosaic protein is Immunoreactive with sera of SARS-infected individuals.

Formulation:

25mM Tris-HCl, 0.4% sarcosyl, 0.25% Triton 100 and 50% glycerol.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

SARS Mosaic antigen is suitable for ELISA and Western blots, excellent antigen for detection of SARS with minimal specificity problems.

Introduction:

SARS Coronavirus is an enveloped virus containing three outer structural proteins, namely the membrane (M), envelope (E), and spike (S) proteins. Spike (S)-glycoprotein of the virus interacts with a cellular receptor and mediates membrane fusion to allow viral entry into susceptible target cells. Accordingly, S-protein plays an important role in virus infection cycle and is the primary target of neutralizing antibodies.

Storage:

SARS Mosaic protein is shipped at ambient temperature. Upon arrival, store at -20C.

To place an order, please [Click HERE](#).